

Anti-Hdac6 (Hd6, Sfc6, Hdac5, Mhda2) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Z2V5

Item Number: CM0012477

Introduction

High-quality HDAC6 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA in mouse and rat. Targets the C-terminal region (aa 836-1037) of mouse HDAC6, a key deacetylase in microtubule dynamics and protein degradation.

[Learn More](#)

Parameter	Value
Product Name	HDAC6 Rabbit mAb
Remarks/Alias	Hd6; Sfc6; Hdac5; mHDA2
Species	Mouse
GeneID (Human)	10013
GeneID	15185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 836-1037 of mouse Hdac6 (Q9Z2V5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Z2V5
Protein Weight	126 kDa
Shipping	Ice bag