


Anti-Hdac6 Monoclonal Antibody For Wb, Elisa - Q9Ubn7

Item Number: CM0007478



Introduction

Anti-HDAC6 rabbit monoclonal antibody validated for WB and ELISA. Recognizes human HDAC6 (~131 kDa), a cytoplasmic deacetylase regulating autophagy, aggresome formation, and microtubule dynamics. Ideal for cancer and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	HDAC6 Rabbit mAb
Remarks/Alias	HD6; JM21; CPBHM; PPP1R90; HDAC6
Species	Human
Gene ID (Human)	10013
Gene ID	10013
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 600-700 of human HDAC6 (NP_006035.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UBN7
Protein Weight	131kDa
Shipping	Ice bag