


Anti-Hdac6 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Ubn7

Item Number: CM0008905



Introduction

Rabbit monoclonal antibody targeting human HDAC6 (Q9UBN7). Validated for WB, IHC-P, IP, ELISA. HDAC6 deacetylase regulates microtubule dynamics, autophagy, and aggresome formation, relevant for cancer and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	HDAC6 Rabbit mAb
Remarks/Alias	HD6; JM21; CPBHM; PPP1R90; HDAC6
Species	Human
GeneID (human)	10013
GeneID	10013
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HDAC6 (Q9UBN7).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	Q9UBN7
Protein Weight	131kDa
Shipping	Ice bag