


Anti-Hdac6 Rabbit Polyclonal Antibody - Q9Ubn7

Item Number: CM0016026



Introduction

Rabbit polyclonal antibody against human HDAC6, validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Suitable for studying tubulin deacetylation, autophagy, and aggresome formation.

[Learn More](#)

Parameter	Value
Product Name	HDAC6 Rabbit pAb
Remarks/Alias	HD6; JM21; CPBHM; PPP1R90; HDAC6
Species	Human
Gene ID (Human)	10013
Gene ID	10013
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 836-1104 of human HDAC6 (NP_006035.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBN7
Protein Weight	131kDa
Shipping	Ice bag