

Anti-Hdac7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Wui4

Item Number: CM0025746



Introduction

Rabbit polyclonal antibody against human HDAC7, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for epigenetic and transcriptional regulation research. - Q8WUI4.

[Learn More](#)

Parameter	Value
Product Name	HDAC7 Rabbit pAb
Remarks/Alias	HD7; HD7A; HDAC7A; HDAC7
Species	Human
Gene ID	51564
Immunogen	Synthetic peptide, A synthetic peptide corresponding to a sequence within amino acids 400-500 of human HDAC7 (XP_011536783.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WUI4
Protein Weight	103kDa
Shipping	Ice bag