


Anti-Hdac8 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9By41

Item Number: CM0029006



Introduction

CamelBio supplies high-quality Anti-HDAC8 Rabbit Polyclonal Antibody (Q9BY41) validated for WB, IHC-P, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, rat. Targets a 42kDa histone deacetylase involved in epigenetic regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HDAC8 Rabbit pAb
Remarks/Alias	HD8; WTS; RPD3; CDA07; CDLS5; KDAC8; MRXS6; HDACL1; C8
Species	Human
Gene ID (Human)	55869
Gene ID	55869
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-377 of human HDAC8 (NP_060956.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BY41
Protein Weight	42kDa
Shipping	Ice bag