


## Anti-Hdac9 Rabbit Polyclonal Antibody For Wb - Q9Ukv0

Item Number: CM0014986



### Introduction

Anti-HDAC9 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. HDAC9 deacetylates histones for epigenetic repression; 111 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | HDAC9 Rabbit pAb                                                                                                         |
| Remarks/Alias    | HD7; HD9; HD7b; HDAC; HDRP; MITR; HDAC7; HDAC7B; HDAC9B; HDAC9FL; HDAC9                                                  |
| Species          | Human                                                                                                                    |
| Gene ID (Human)  | 9734                                                                                                                     |
| Gene ID          | 9734                                                                                                                     |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HDAC9 (NP_478056.1). |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, IHC-P, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | Q9UKV0                                                                                                                   |
| Protein Weight   | 111kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                  |