

Anti-Hdc Rabbit Polyclonal Antibody For Wb, Elisa - P19113

Item Number: CM0029122



Introduction

High-quality rabbit polyclonal antibody targeting Histidine decarboxylase (HDC), validated for Western blot (WB) and ELISA applications. Reacts with human and mouse samples. Ideal for studying histamine biosynthesis and catecholamine metabolism.

[Learn More](#)

Parameter	Value
Product Name	HDC Rabbit pAb
Remarks/Alias	HDC
Species	Human
Gene ID	3067 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 453-662 of human HDC (NP_002103.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19113
Protein Weight	74kDa
Shipping	Ice bag