


Anti-Hdgf Monoclonal Antibody For Wb, If-P, Elisa - P51858

Item Number: CM0002958



Introduction

Rabbit monoclonal HDGF antibody for Western blot, immunofluorescence, and ELISA. Validated in human, mouse, and rat. Recognizes hepatoma-derived growth factor (HMG1L2), a transcriptional repressor with mitogenic and heparin-binding activities—suitable for cancer and growth factor research.

[Learn More](#)

Parameter	Value
Product Name	HDGF Rabbit mAb
Remarks/Alias	HMG1L2; HDGF
Species	Human
GeneID (Human)	3068
GeneID	3068
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HDGF (P51858).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51858
Protein Weight	27kDa
Shipping	Ice bag