


Anti-Hdgf Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P51858

Item Number: CM0028776



Introduction

Rabbit polyclonal antibody targeting HDGF (HMG1L2), validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Suitable for studying transcriptional repression, mitogenic activity, and heparin-binding functions.

[Learn More](#)

Parameter	Value
Product Name	HDGF Rabbit pAb
Remarks/Alias	HMG1L2; HDGF
Species	Human
GeneID (Human)	3068
GeneID	3068
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human HDGF (NP_004485.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51858
Protein Weight	27kDa
Shipping	Ice bag