


Anti-Hdlbp Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q00341

Item Number: CM0029794



Introduction

HDLBP (Vigilin) rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat HDLBP (141 kDa). Ideal for cholesterol metabolism and lipoprotein research.

[Learn More](#)

Parameter	Value
Product Name	HDLBP Rabbit pAb
Remarks/Alias	HBP; VGL; PRO2900; HDLBP
Species	Human
GeneID (human)	3069
GeneID	3069
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1035-1268 of human HDLBP (NP_976221.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00341
Protein Weight	141kDa
Shipping	Ice bag