


Anti-Hdlbp/Vigilin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q00341

Item Number: CM0003738



Introduction

Rabbit monoclonal antibody against human HDLBP (Vigilin), validated in WB, IHC-P, and ELISA. Detects endogenous levels of HDLBP in human, mouse, and rat samples. Ideal for lipid metabolism and cholesterol research.

[Learn More](#)

Parameter	Value
Product Name	HDLBP Rabbit mAb
Remarks/Alias	HBP; VGL; PRO2900; HDLBP
Species	Human
GeneID (human)	3069
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HDLBP (NP_976221.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00341
Protein Weight	141kDa
Shipping	Ice bag