

Anti-He4/Wfdc2 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q14508

Item Number: CM0008769

Introduction

High-quality rabbit monoclonal antibody targeting human HE4/WFDC2 (Q14508). Validated for WB, IF/ICC, ELISA and cross-reacts with mouse. Ideal for studying lung development, reproductive biology, and cancer biomarkers.

[Learn More](#)

Parameter	Value
Product Name	HE4/WFDC2 Rabbit mAb
Remarks/Alias	HE4; WAP5; EDDM4; dJ461P17.6; HE4/WFDC2
Species	Human
GeneID (Human)	10406
GeneID	10406
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 45-124 of human HE4/WFDC2 (Q14508).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14508
Protein Weight	13kDa
Shipping	Ice bag