


Anti-Hec1/Ndc80 Polyclonal Antibody For Wb, Elisa - O14777

Item Number: CM0030077



Introduction

Rabbit polyclonal antibody targeting human HEC1/NDC80, validated for western blot and ELISA. Suitable for kinetochore and mitosis research, with cross-reactivity to mouse.

[Learn More](#)

Parameter	Value
Product Name	HEC1/NDC80 Rabbit pAb
Remarks/Alias	HEC; HEC1; TID3; KNTC2; HsHec1; hsNDC80; HEC1/NDC80
Species	Human
GeneID (Human)	10403
GeneID	10403
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 343-642 of human HEC1/NDC80 (NP_006092.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14777
Protein Weight	74kDa
Shipping	Ice bag