


Anti-Hec1/Ndc80 Rabbit Monoclonal Antibody For Wb, Elisa - O14777

Item Number: CM0007112



Introduction

Rabbit monoclonal antibody against human HEC1/NDC80 (O14777). Validated for Western blot (WB) and ELISA. Ideal for kinetochore function, mitotic checkpoint, and cancer studies. Predicted molecular weight: 74 kDa.

[Learn More](#)

Parameter	Value
Product Name	HEC1/NDC80 Rabbit mAb
Remarks/Alias	HEC; HEC1; TID3; KNTC2; HsHec1; hsNDC80; HEC1/NDC80
Species	Human
GeneID	10403
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 543-642 of human HEC1/NDC80 (O14777)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14777
Protein Weight	74kDa
Shipping	Ice bag