

Anti-Hectd1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ult8

Item Number: CM0024719



Introduction

Rabbit polyclonal antibody targeting human HECTD1 (EULIR), suitable for WB and ELISA. Detects endogenous HECTD1 at ~289kDa in human samples. Ideal for research on ubiquitination, neural development, and stem cell biology.

[Learn More](#)

Parameter	Value
Product Name	HECTD1 Rabbit pAb
Remarks/Alias	EULIR; HECTD1
Species	Human
GeneID (Human)	25831
GeneID	25831
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2518-2610 of human HECTD1 (NP_056197.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9ULT8
Protein Weight	289kDa
Shipping	Ice bag