

Anti-Hel308 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Tdg4

Item Number: CM0015364

Introduction

Rabbit polyclonal antibody against human HEL308 helicase for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for DNA double-strand break repair, microhomology-mediated end-joining, and RAD51 pathway studies.

[Learn More](#)

Parameter	Value
Product Name	HEL308 Rabbit pAb
Remarks/Alias	HEL308
Species	Human
Gene ID (Human)	113510
Gene ID	113510
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human HEL308 (NP_598375.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDG4
Protein Weight	124kDa
Shipping	Ice bag