

Anti-Heme Oxygenase 1 (Ho-1/Hmox1) Rabbit Polyclonal Antibody - P09601

Item Number: CM0014843



Introduction

Rabbit polyclonal antibody against human Heme Oxygenase 1 (HO-1/HMOX1). Validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Suitable for apoptosis and heme catabolism research.

[Learn More](#)

Parameter	Specification
Product Name	Heme Oxygenase 1 (HO-1/HMOX1) Rabbit pAb
Remarks/Alias	HO-1; HSP32; HMOX1D; bK286B10; Heme Oxygenase 1 (HO-1/HMOX1)
Species	Human
GeneID	3162
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-288 of human Heme Oxygenase 1 (HO-1/HMOX1) (NP_002124.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09601
Protein Weight	33kDa
Shipping	Ice bag