

Anti-Heme Oxygenase 1 (Ho-1/Hmox1) Rabbit Polymab For Wb, Ip, If, Ihc, Elisa - P09601 / P14901

Item Number: CM0012660



Introduction

Rabbit recombinant monoclonal antibody targeting Heme Oxygenase 1 (HO-1/HMOX1). Validated for WB, IP, IF, IHC, and ELISA applications in human, mouse, and rat samples. Detects a 33 kDa protein. Suitable for apoptosis, heme metabolism, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Heme Oxygenase 1 (HO-1/HMOX1) Rabbit PolymAb®
Remarks/Alias	HO-1; HSP32; HMOX1D; bK286B10; D8Wsu38e; HO-1; HO1; Hemox; Hmox; Hsp32
Species	Human/Mouse
Gene ID (Human)	3162
Gene ID	3162/15368
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1-266 of mouse Hmox1 (NP_034572.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09601/P14901
Protein Weight	33 kDa
Shipping	Ice bag