

Anti-Heme Oxygenase 1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P14901

Item Number: CM0007819

Introduction

Rabbit monoclonal antibody targeting Heme Oxygenase 1 (HO-1/Hmox1). Validated for WB, IF-P, IHC-P, and ELISA applications. Detects endogenous HO-1 in human, mouse, and rat samples. Useful for oxidative stress and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	Heme Oxygenase 1 Rabbit mAb
Remarks/Alias	HO1; HO-1; Hmox; Hemox; Hsp32; D8Wsu38e
Species	Mouse
Gene ID (Human)	3162
Gene ID	15368
Immunogen	Recombinant protein corresponding to amino acids 1-266 of mouse Hmox1 (NP_034572.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14901
Protein Weight	33 kDa
Shipping	Ice bag