

Anti-Hemgn Rabbit Polyclonal Antibody For Wb, Elisa - Q9BxL5

Item Number: CM0021126

Introduction

Rabbit polyclonal antibody targeting human HEMGN (Hemogen) for WB and ELISA applications. Cross-reacts with human and mouse; suitable for hematopoietic differentiation, leukemia, and NF-κB apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	HEMGN Rabbit pAb
Remarks/Alias	NDR; EDAG; CT155; EDAG-1; HEMGN
Species	Human
GeneID (Human)	55363
GeneID	55363
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human HEMGN (NP_932095.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BXL5
Protein Weight	55kDa
Shipping	Ice bag