


Anti-Hemoglobin Subunit Alpha (Hba1) Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P69905

Item Number: CM0007887



Introduction

High-quality rabbit monoclonal antibody targeting human Hemoglobin subunit alpha (HBA1). Validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for oxygen transport and hemoglobinopathies research.

[Learn More](#)

Parameter	Value
Product Name	Hemoglobin subunit alpha (HBA1) Rabbit mAb
Remarks/Alias	HBH; ECT7; HBA-T3; METHBA; Hemoglobin subunit alpha (HBA1)
Species	Human
GeneID	3039
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HBA1 (P69905)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P69905
Protein Weight	15kDa
Shipping	Ice bag