


Anti-Hemoglobin Subunit Alpha (Hba1) Rabbit Polyclonal Antibody For Wb And Elisa - P69905

Item Number: CM0016786



Introduction

Rabbit polyclonal antibody targeting human Hemoglobin subunit alpha (HBA1) for Western blot and ELISA. Ideal for oxygen transport, hemoglobinopathy, and red blood cell research.

[Learn More](#)

Parameter	Value
Product Name	Hemoglobin subunit alpha (HBA1) Rabbit pAb
Remarks/Alias	HBH; ECYT7; HBA-T3; METHBA; Hemoglobin subunit alpha (HBA1)
Species	Human
GeneID (Human)	3039
GeneID	3039
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Hemoglobin subunit alpha (NP_000549.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P69905
Protein Weight	15kDa
Shipping	Ice bag