

## Anti-Hemopexin (Hpx) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P02790

Item Number: CM0008634

### Introduction

Rabbit monoclonal antibody against human Hemopexin (HPX), validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with rat. Essential for heme transport and iron recovery studies. Ideal for diagnostic and research applications.

[Learn More](#)

| Parameter        | Value                                                                                                        |
|------------------|--------------------------------------------------------------------------------------------------------------|
| Product Name     | Hemopexin (HPX) Rabbit mAb                                                                                   |
| Remarks/Alias    | HX; Hemopexin (HPX)                                                                                          |
| Species          | Human                                                                                                        |
| GeneID (human)   | 3263                                                                                                         |
| GeneID           | 3263                                                                                                         |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 363-462 of human Hemopexin (HPX) (P02790) |
| Source           | Rabbit                                                                                                       |
| Category         | Monoclonal Antibodies                                                                                        |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                       |
| Cross Reactivity | Human, Rat                                                                                                   |
| SWISS            | P02790                                                                                                       |
| Protein Weight   | 52kDa                                                                                                        |
| Shipping         | Ice bag                                                                                                      |