

Anti-Hepacam Rabbit Polyclonal Antibody For Wb, Elisa - Q14Cz8

Item Number: CM0016641



Introduction

Rabbit polyclonal antibody targeting human HEPACAM (46 kDa), validated for WB and ELISA. Reacts with human, mouse, and rat species. Ideal for investigating cell adhesion, migration, and growth regulation.

[Learn More](#)

Parameter	Value
Product Name	HEPACAM Rabbit pAb
Remarks/Alias	HEPN1; MLC2A; MLC2B; GlialCAM; HEPACAM
Species	Human
Gene ID (Human)	220296
Gene ID	220296
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-245 of human HEPACAM (NP_689935.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14CZ8
Protein Weight	46kDa
Shipping	Ice bag