

Anti-Heparanase 1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y251

Item Number: CM0020080



Introduction

Rabbit polyclonal antibody against human Heparanase 1 (HPSE), suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for ECM degradation, tumor invasion, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Heparanase 1 Rabbit pAb
Remarks/Alias	HPA; HPA1; HPR1; HSE1; HPSE1; Heparanase 1
Species	Human
Gene ID (Human)	10855
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-275 of human Heparanase 1 (NP_001092010.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y251
Protein Weight	61kDa
Shipping	Ice bag