


Anti-Her2/ErbB2 Rabbit Monoclonal Antibody For Wb, Ihc, If - P04626

Item Number: CM0002338



Introduction

Recombinant rabbit monoclonal antibody against human HER2/ErbB2 (P04626), validated for WB, IHC, IF, and ELISA. Ideal for cancer signaling research and diagnostic assay development.

[Learn More](#)

Parameter	Specification
Product Name	HER2/ErbB2 Rabbit mAb
Remarks/Alias	NEU; NGL; HER2; TKR1; CD340; HER-2; VSCN2; MLN 19; MLN-19; c-ERB2; c-ERB-2; HER-2/neu; p185(erbB2)
Species	Human
GeneID	2064
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-652 of human HER2/ErbB2 (P04626)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, mIHC, ELISA
Cross Reactivity	Human
SWISS	P04626
Protein Weight	138kDa
Shipping	Ice bag