

Anti-Her2/ErbB2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04626

Item Number: CM0021800



Introduction

Rabbit polyclonal antibody against HER2/ErbB2 (receptor tyrosine kinase), validated for WB, IHC-P, ELISA, reacts with human, mouse, rat. Suitable for cancer research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	HER2/ErbB2 Rabbit pAb
Remarks/Alias	NEU; NGL; HER2; TKR1; CD340; HER-2; VSCN2; MLN 19; MLN-19; c-ERB2; c-ERB-2; HER-2/neu; p185(erbB2); HER2/ErbB2
Species	Human
GeneID (Human)	2064
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1156-1255 of human HER2/ErbB2 (NP_004439.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04626
Protein Weight	138kDa
Shipping	Ice bag