

Anti-Herc3 Rabbit Polyclonal Antibody For Wb, Elisa - Q15034

Item Number: CM0030684

Introduction

HERC3 rabbit polyclonal antibody for WB and ELISA; reacts with human and mouse.

Immunogen: Recombinant fusion protein corresponding to amino acids 1-200 of human HERC3 (NP_055421.1).

[Learn More](#)

Parameter	Value
Product Name	HERC3 Rabbit pAb
Remarks/Alias	HERC3
Species	Human
Gene ID (Human)	8916
Gene ID	8916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human HERC3 (NP_055421.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15034
Protein Weight	117kDa
Shipping	Ice bag