


Anti-Hes1 Polyclonal Antibody For Wb And Elisa - Q14469

Item Number: CM0013022



Introduction

High-quality rabbit anti-HES1 polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat HES1, a transcriptional repressor involved in Notch signaling and development. Ideal for gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	HES1 Rabbit pAb
Remarks/Alias	HHL; HRY; HES-1; bHLHb39; HES1
Species	Human
Gene ID (Human)	3280
Gene ID	3280
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 90-165 of human HES1 (NP_005515.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14469
Protein Weight	30kDa
Shipping	Ice bag