


Anti-Hes1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q14469

Item Number: CM0003263



Introduction

Rabbit monoclonal antibody targeting human HES1 transcription factor. Suitable for Western blot, IHC-P, IP, and ELISA. Reacts with human, mouse, and rat. Research use only.

[Learn More](#)

Parameter	Value
Product Name	HES1 Rabbit mAb
Remarks/Alias	HHL; HRY; HES-1; bHLHb39; HES1
Species	Human
GeneID	3280
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 181-280 of human HES1 (Q14469)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14469
Protein Weight	30kDa
Shipping	Ice bag