


## Anti-Hes1 Rabbit Polyclonal Antibody For Wb, Elisa - Q14469

Item Number: CM0017461



### Introduction

Rabbit polyclonal antibody against human transcription factor HES1, validated for WB and ELISA. Detects human HES1 with cross-reactivity to rat and mouse. Suitable for transcriptional regulation and developmental biology research.

[Learn More](#)

| Parameter        | Value                                                                                                                                |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | HES1 Rabbit pAb                                                                                                                      |
| Remarks/Alias    | HHL; HRY; HES-1; bHLHb39; HES1                                                                                                       |
| Species          | Human                                                                                                                                |
| GeneID (Human)   | 3280                                                                                                                                 |
| GeneID           | 3280                                                                                                                                 |
| Immunogen        | synthetic peptide Recombinant fusion protein containing a sequence corresponding to amino acids 151-280 of human HES1 (NP_005515.1). |
| Source           | Rabbit                                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                                |
| Application      | WB, ELISA                                                                                                                            |
| Cross Reactivity | Rat, Mouse                                                                                                                           |
| SWISS            | Q14469                                                                                                                               |
| Protein Weight   | 30kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                              |