


Anti-Hex Rabbit Monoclonal Antibody For Wb, Elisa - Q03014

Item Number: CM0002950



Introduction

Rabbit monoclonal antibody against human HHEX/Hex (Q03014) validated for WB and ELISA. A transcriptional regulator in WNT signaling, hematopoietic differentiation, and anterior development. Ideal for transcription studies and assay development.

[Learn More](#)

Parameter	Value
Product Name	Hex Rabbit mAb
Remarks/Alias	HEX; PRH; HMPH; PRHX; HOX11L-PEN; Hex
Species	Human
GeneID (human)	3087
GeneID	3087
Immunogen	Synthetic peptide corresponding to amino acids 50-150 of human Hex (Q03014)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q03014
Protein Weight	30kDa
Shipping	Ice bag