


Anti-Hexa (Beta-Hexosaminidase Subunit Alpha) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P06865

Item Number: CM0006588



Introduction

Rabbit monoclonal antibody against human HEXA, validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying lysosomal storage disorders like Tay-Sachs disease and glycosphingolipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	HEXA Rabbit mAb
Remarks/Alias	TSD; HEXA
Species	Human
GeneID (human)	3073
GeneID	3073
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 368-529 of human HEXA. (NP_000511.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06865
Protein Weight	61kDa
Shipping	Ice bag