


Anti-Hexim1 Rabbit Monoclonal Antibody - O94992

Item Number: CM0006487



Introduction

Anti-HEXIM1 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, IP, and ELISA. Detects human, mouse, and rat HEXIM1 (~41 kDa), a transcriptional regulator and P-TEFb inhibitor.

[Learn More](#)

Parameter	Specification
Product Name	HEXIM1 Rabbit mAb
Remarks/Alias	CLP1, EDG1, HIS1, MAQ1, HEXIM1
Species	Human
Gene ID (Human)	10614
Gene ID	10614
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 255-359 of human HEXIM1 (NP_006451.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94992
Protein Weight	41kDa
Shipping	Ice bag