


Anti-Hexim1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 094992

Item Number: CM0028492



Introduction

High-quality rabbit polyclonal antibody targeting human HEXIM1, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for transcription regulation and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	HEXIM1 Rabbit pAb
Remarks/Alias	CLP1; EDG1; HIS1; MAQ1; HEXIM1
Species	Human
Gene ID (Human)	10614
Gene ID	10614
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-359 of human HEXIM1 (NP_006451.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	094992
Protein Weight	41kDa
Shipping	Ice bag