

Anti-Hexim2 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Mh2

Item Number: CM0018364

Introduction

Rabbit polyclonal antibody targeting human HEXIM2 (UniProt Q96MH2), validated in WB and ELISA. Cross-reacts with mouse. Essential tool for transcription regulation and 7SK RNP complex studies.

[Learn More](#)

Parameter	Value
Product Name	HEXIM2 Rabbit pAb
Remarks/Alias	L3; HEXIM2
Species	Human
Gene ID (Human)	124790
Gene ID	124790
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 147-286 of human HEXIM2 (NP_653209.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96MH2
Protein Weight	32kDa
Shipping	Ice bag