

Anti-Hexokinase II Polyclonal Antibody For Wb, Elisa - P52789

Item Number: CM0013415



Introduction

Rabbit polyclonal antibody against human Hexokinase II, validated for WB and ELISA in human, mouse, and rat. Targets the glycolytic enzyme HK2, suitable for metabolism and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Hexokinase II Rabbit pAb
Remarks/Alias	HKII, HXK2, II
Species	Human
GeneID (Human)	3099
GeneID	3099
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1-120 of human Hexokinase II (NP_000180.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52789
Protein Weight	102kDa
Shipping	Ice bag