

Anti-Hexokinase II Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P52789

Item Number: CM0014978



Introduction

Rabbit polyclonal antibody against human Hexokinase II (HKII/HXK2) validated for Western blot, immunohistochemistry, immunofluorescence, and ELISA; detects human, mouse, and rat. Suitable for glycolysis and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Hexokinase II Rabbit pAb
Remarks/Alias	HKII; HXK2; Hexokinase II
Species	Human
Gene ID	3099
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human Hexokinase II (NP_000180.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52789
Protein Weight	102kDa
Shipping	Ice bag