


Anti-Hey1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y5J3

Item Number: CM0024153



Introduction

High-quality rabbit polyclonal antibody against human HEY1, validated for Western blot, IF/ICC, and ELISA. Ideal for cardiovascular development, Notch signaling, and transcriptional repression research.

[Learn More](#)

Parameter	Value
Product Name	HEY1 Rabbit pAb
Remarks/Alias	CHF2; OAF1; HERP2; HESR1; HRT-1; NERP2; hHRT1; BHLHb31; HEY1
Species	Human
Gene ID (Human)	23462
Gene ID	23462
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 114-129 of human HEY1 (NP_001269780.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5J3
Protein Weight	33kDa
Shipping	Ice bag