

Anti-Hey2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ubp5

Item Number: CM0012867

Introduction

High-quality HEY2 rabbit polyclonal antibody (pAb) for western blot and ELISA, reactive against human, mouse, rat. Suitable for studying cardiovascular development, Notch signaling, and transcriptional repression mechanisms.

[Learn More](#)

Parameter	Value
Product Name	HEY2 Rabbit pAb
Remarks/Alias	GRL; CHF1; HRT2; HERP1; HESR2; bHLHb32; GRIDLOCK; HEY2
Species	Human
GeneID (Human)	23493
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human HEY2 (NP_036391.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBP5
Protein Weight	36kDa
Shipping	Ice bag