

Anti-Heyl Rabbit Polyclonal Antibody (Hey3, Hrt3, Bhlhb33) For Wb, Elisa - Q9Nq87

Item Number: CM0024026

Introduction

HEYL Rabbit Polyclonal Antibody (pAb) validated for WB and ELISA in human and mouse. Targets HEYL transcriptional repressor in Notch signaling, key for cardiovascular development and gene regulation. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	HEYL Rabbit pAb
Remarks/Alias	HEY3; HRT3; HESR3; bHLHb33; HEYL
Species	Human
GeneID (Human)	26508
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human HEYL (NP_055386.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NQ87
Protein Weight	35kDa
Shipping	Ice bag