

Anti-Hfe Monoclonal Antibody For Wb, Elisa - Q30201

Item Number: CM0002416



Introduction

This rabbit monoclonal antibody targets human HFE (Hereditary Hemochromatosis Protein), validated for Western blotting and ELISA applications. Suitable for investigating iron metabolism, hemochromatosis, and transferrin receptor interactions.

[Learn More](#)

Parameter	Value
Product Name	HFE Rabbit mAb
Remarks/Alias	HH; HFE1; HLA-H; MVCD7; TFQTL2; HFE
Species	Human
GeneID (human)	3077
GeneID	3077
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 249-348 of human HFE (Q30201).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q30201
Protein Weight	40kDa
Shipping	Ice bag