

Anti-Hfe Polyclonal Antibody For Wb, Elisa - Q30201

Item Number: CM0017424



Introduction

High-quality rabbit polyclonal antibody targeting human HFE (hereditary hemochromatosis protein), validated for Western blot and ELISA. Reacts with human and mouse samples. Supports iron metabolism research.

[Learn More](#)

Parameter	Value
Product Name	HFE Rabbit pAb
Remarks/Alias	HH; HFE1; HLA-H; MVCD7; TFQTL2; HFE
Species	Human
Gene ID (Human)	3077
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-306 of human HFE (NP_000401.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q30201
Protein Weight	40kDa
Shipping	Ice bag