

Anti-Hfm1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - A2Pyh4

Item Number: CM0025137

Introduction

This rabbit polyclonal antibody specifically recognizes human HFM1 and cross-reacts with mouse and rat, validated for WB, IHC-P, and ELISA. Ideal for meiosis and premature ovarian failure studies.

[Learn More](#)

Parameter	Value
Product Name	HFM1 Rabbit pAb
Remarks/Alias	MER3; POF9; Si-11; SEC63D1; Si-11-6; helicase; HFM1
Species	Human
GeneID	164045
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human HFM1 (NP_001017975.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	A2PYH4
Protein Weight	163kDa
Shipping	Ice bag