


Anti-Hgd Rabbit Monoclonal Antibody For Wb, Elisa - Q93099

Item Number: CM0003921



Introduction

Rabbit monoclonal antibody targeting human HGD (Homogentisate 1,2-dioxygenase) for WB and ELISA. Ideal for metabolic disease research, including alkaptonuria. Specificity ensured by peptide immunogen.

[Learn More](#)

Parameter	Value
Product Name	HGD Rabbit mAb
Remarks/Alias	AKU; HGO; HGD
Species	Human
GeneID (Human)	3081
GeneID	3081
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 346-445 of human HGD (Q93099).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q93099
Protein Weight	50kDa
Shipping	Ice bag