


Anti-Hgfac Polyclonal Antibody For Wb And Elisa - Q04756

Item Number: CM0018259



Introduction

HGFAC Rabbit pAb for WB and ELISA. Rabbit polyclonal antibody against human HGFAC (aa 408-655). Cross-reacts with mouse and rat. Ideal for HGF activation studies.

[Learn More](#)

Parameter	Value
Product Name	HGFAC Rabbit pAb
Remarks/Alias	HGFA; HGFAC
Species	Human
Gene ID (Human)	3083
Gene ID	3083
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 408-655 of human HGFAC (NP_001519.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04756
Protein Weight	71kDa
Shipping	Ice bag