


Anti-Hgs Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - O14964

Item Number: CM0024578



Introduction

Rabbit polyclonal antibody against human HGS (Hepatocyte growth factor-regulated tyrosine kinase substrate), validated for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for studying endosomal sorting and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	HGS Rabbit pAb
Remarks/Alias	HRS; HGS
Species	Human
Gene ID	9146 (Human)
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14964
Protein Weight	86kDa
Shipping	Ice bag