

Anti-Hhla2 Rabbit Polyclonal Antibody For Wb, Ihc-P And Elisa - Q9Um44

Item Number: CM0013523

Introduction

High-quality anti-HHLA2 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA applications. Detects human HHLA2, a T-cell costimulatory ligand, ideal for immuno-oncology research and diagnostic development.

[Learn More](#)

Parameter	Details
Product Name	HHLA2 Rabbit pAb
Remarks/Alias	B7y; B7H7; B7-H5; B7-H7; HHLA2
Species	Human
GeneID (Human)	11148
GeneID	11148
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-200 of human HHLA2 (NP_009003.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9UM44
Protein Weight	47kDa
Shipping	Ice bag