


Anti-Hibadh Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P31937

Item Number: CM0007630



Introduction

HIBADH Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat 3-hydroxyisobutyrate dehydrogenase, a mitochondrial enzyme in BCAA catabolism. Ideal for metabolic research.

[Learn More](#)

Parameter	Specification
Product Name	HIBADH Rabbit mAb
Remarks/Alias	NS5ATP1; HIBADH
Species	Human
GeneID (Human)	11112
GeneID	11112
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-336 of human HIBADH (NP_689953.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31937
Protein Weight	35kDa
Shipping	Ice bag